



**SUN BRAND**

Securely *Sealing* Sentiments  
An ISO CERTIFIED COMPANY

# Adhesive Specialities

*Manufacturer of Speciality Adhesive Tapes*

TECHNICAL DATA SHEET

**Product Code – AS-033-A**

**DOUBLE SIDE TISSUE TAPE – 60 mic – HOTMELT**

**Description:** Double Side Tissue Tape is coated with hotmelt adhesive on both the sides and comes with an easy-to-peel liner.

**Application:** Double Side Tissue Tape is used for sticking poly bag and paper, for holding and joining light weight objects, fixing of ornament hook or signage board, shoe & leather industry and other general applications.

**Application Technique:** For best results, ensure the surface is clean, dry and dust free before application of the tape. Apply even and firm pressure for better adhesion results.

## TECHNICAL VALUES

PROPERTIES	VALUES
Base Material	Tissue
Adhesive	Hotmelt
Total Thickness	60 ± 5 µm
Initial Tack (mm) (Steel Ball #14)	120 mm
Peel Adhesion (g/25mm)	800 gms
Temperature Resistance (°C/hr.)	Adhesive: 60°C
Holding Power (hr.)	> 2
Colour	Translucent
Release Liner	Paper

**Storage & Shelf Life:** The tape should be stored in a cool and dry place and prolonged exposure to the direct sunlight should be avoided. To obtain best performance, the product should be used within 12 months from the manufacturing date.

[www.adhesivespecialities.co.in](http://www.adhesivespecialities.co.in)

©Adhesive Specialities

The technical values are determined by/based on standard tests. Our recommendations on the use of the product are based on reliable in-house tests, but the users are requested to conduct their own test to determine the suitability of the product. Adhesive Specialities will not be responsible for any loss or damage caused as a result of our recommendations.